



Frequently Asked Questions

- **How Does AcceleDent Aura Work?**

By inserting the Mouthpiece – fitted around your existing orthodontics – and turning on the Activator for 20 minutes every day, AcceleDent Aura can speed up tooth movement by using SoftPulse Technology™. Faster tooth movement may decrease the duration of your orthodontic treatment – and can even help to make your orthodontics more comfortable.

- **Is AcceleDent Aura safe?**

Yes. AcceleDent with SoftPulse Technology™ has been demonstrated safe and reliable in US clinical trials.

- **How is the orthodontic treatment accelerated?**

Your orthodontics work on their own by moving your teeth in specific directions. AcceleDent Aura uses SoftPulse Technology™ to generate small vibrations – micropulses – enhancing the movements directed by your orthodontics allowing for faster tooth movement.

- **How do I know if I am biting into the Mouthpiece correctly and with the right amount of force?**

Bite pressure on the Mouthpiece should be just firm enough to hold AcceleDent Aura in place without using your hands. Avoid biting forcibly on the Mouthpiece as that may cause excessive wear.

- **Can I try to wear AcceleDent Aura while I am sleeping?**

No. It is not recommended that AcceleDent Aura be used while sleeping.

- **Can I wear my AcceleDent Aura when it is charging or connected to a computer?**

No, do not wear your AcceleDent Aura when it is charging or connected to a computer.

- **Do I have to do the full 20 minute treatment at one time?**

Yes. It is recommended that you use your AcceleDent Aura for a continuous 20 minute period each day.

- **If I absolutely have to interrupt the 20 minute session for a few minutes, when I turn my AcceleDent Aura back on, does it remember how long I used it?**

Yes. AcceleDent Aura will remember how long you used it as long as the interruption does not exceed 5 minutes. After 5 minutes of being turned off, AcceleDent Aura will work for a full session of 20 minutes at its next use.

- **If I use AcceleDent Aura for more than the recommended 20 minute treatment time, does it add incremental benefit?**

No. Additional daily use greater than the prescribed 20 minutes duration is not recommended.

- **How often does AcceleDent Aura need to be recharged?**

The Activator should be recharged every 3 days.

- **What is the best way to clean AcceleDent Aura?**

The Mouthpiece should first be disconnected from the Activator and then rinsed well with lukewarm water after each use. You may occasionally use a soft liquid dishwashing detergent, a mild soap, or toothpaste to gently brush the Mouthpiece. The Activator can be wiped with a soft cloth as necessary. It is

recommended that you allow AcceleDent to air dry before being stored in the Travel Case. **Do not place the Mouthpiece or the Activator in a dishwasher as this will void the warranty and could cause damage to the product.**

- **Can I use my AcceleDent Aura in the shower?**

No. AcceleDent Aura is not water-resistant and should not be used in the shower.

- **Can I sanitize my AcceleDent Aura and share it with a friend or family member who also has braces?**

No. AcceleDent Aura is an FDA cleared, single-user prescription device. Sharing the AcceleDent Aura and/or the Mouthpiece is prohibited. Your orthodontist prescribed AcceleDent Aura for your use only, based upon your unique orthodontic treatment plan. Use of AcceleDent Aura without the supervision and direction of an orthodontist may result in adverse unintended consequences.

- **What type of battery does AcceleDent Aura contain?**

The rechargeable battery in AcceleDent Aura is a medical grade Lithium ion battery.

- **What shall I do with my AcceleDent Aura once treatment is finished?**

You should dispose of your AcceleDent Aura as you would with any other device with a rechargeable battery. Reference the “Disposal” section of the Directions for Use manual for further details.

- **Is my AcceleDent Aura under warranty?**

Yes, the details of the warranty can be found in the DFU manual, click [here](#).

- **What does it mean if my Activator stops working and the LED Indicator is flashing alternating green and orange?**

The Activator will deactivate after completing the equivalent of 18 months of 20 minute daily sessions. If you are in the middle of a treatment plan please contact your orthodontist.

- **Why did the date change on my FastTrac Usage Report?**

If your Activator battery is fully drained and the activator shuts down (no flashing LED Indicator or vibrations), the date on the Usage report will reset to January 1, 2012. The battery has to be charged when the LED indicates a low battery (shown when the LED Indicator is quickly flashing orange and the Activator is powered on but not charging). Refer to the FastTrac Usage Report section in the DFU to reset date, click [here](#).