

CASE STUDY



DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS

INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE

TREATMENT TIME

PROJECTED TREATMENT TIME:
13 MONTHS

ACTUAL TREATMENT TIME WITH
ACCELEDDENT:
5 MONTHS

ADDITIONAL RESOURCES

ON DEMAND CE PROGRAM

[Accelerated Aligner Therapy: Precision Taken to the Next Level](#)

Presenter: Dr. Manal Ibrahim

CASE STUDY



DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS

INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
13 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
5 MONTHS



BASELINE



CASE STUDY



CLINCHECK® SOFTWARE TREATMENT INFORMATION

DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS


INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
13 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
5 MONTHS


TOOTH MOVEMENT ASSESMENT



TOOTH MOVEMENT ASSESSMENT

Tooth	Assessment	Movement	Value	Range
UL4	Blue	Extrusion	0.9 mm	0.5 mm-1.0 mm
UL5	Black	Extrusion	1.1 mm	> 1.0 mm
LL7	Blue	Extrusion	0.7 mm	0.5 mm-1.0 mm
LR2	Blue	Intrusion	2.6 mm	2.5 mm-3.0 mm
LR3	Blue	Root Movement	4.7 mm	4.0 mm-6.0 mm
LR4	Blue	Intrusion	0.6 mm	0.5 mm-1.0 mm

STAGING



CASE STUDY



DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS

INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
13 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
5 MONTHS

PROGRESS

3 Months



CASE STUDY



DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS

INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
13 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
5 MONTHS



FINISH



CASE STUDY



DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS

INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
13 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
5 MONTHS

BASELINE





PROGRESS 3 MONTHS



FINISH 5 MONTHS



CLINCHECK SOFTWARE PROJECTED OUTCOME



CASE STUDY



CLINCHECK
SOFTWARE
PROJECTED
OUTCOME

DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS

INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
13 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
5 MONTHS

BASELINE



PROGRESS 3 MONTHS



FINISH 5 MONTHS



CASE STUDY



DEMOGRAPHICS

FEMALE
27 YEARS OLD

DIAGNOSIS

CLASS I
SPACING

MECHANICS

INVISALGN®
26 TOTAL ALIGNERS
NO REFINEMENT
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
13 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
5 MONTHS

