

# CASE STUDY



DEMOGRAPHICS	DIAGNOSIS	MECHANICS	TREATMENT TIME
MALE 14 YEARS OLD	CLASS I, II DIV 2 CROWDING DEEP BITE (50%) NARROW ARCHES EXCESSIVE OVERJET (5MM)	INVISALIGN® 33 TOTAL ALIGNERS NO REFINEMENT 7 DAY ALIGNER CHANGE	PROJECTED TREATMENT TIME: 18 MONTHS  ACTUAL TREATMENT TIME WITH ACCELEDENT: 8 MONTHS

*“Another doctor had recommended headgear & braces for 2 years. We offered Invisalign & AcceleDent and finished in 8 months.”*

- Dr. Jasmine Gorton

# CASE STUDY



## DEMOGRAPHICS

MALE  
14 YEARS OLD

## DIAGNOSIS

CLASS I, II DIV 2  
CROWDING  
DEEP BITE (50%)  
NARROW ARCHES  
EXCESSIVE OVERJET (5MM)

## MECHANICS

INVISALIGN®  
33 TOTAL ALIGNERS  
NO REFINEMENT  
7 DAY ALIGNER CHANGE

## TREATMENT TIME

PROJECTED TREATMENT  
TIME:  
18 MONTHS

ACTUAL TREATMENT  
TIME WITH ACCELEDENT:  
8 MONTHS



## BASELINE

“Another doctor had recommended headgear & braces for 2 years. We offered Invisalign & AcceleDent and finished in 8 months.”



# CASE STUDY



## DEMOGRAPHICS

MALE  
14 YEARS OLD

## DIAGNOSIS

CLASS I, II DIV 2  
CROWDING  
DEEP BITE (50%)  
NARROW ARCHES  
EXCESSIVE OVERJET (5MM)

## MECHANICS

INVISALIGN®  
33 TOTAL ALIGNERS  
NO REFINEMENT  
7 DAY ALIGNER CHANGE

## TREATMENT TIME

PROJECTED TREATMENT  
TIME:  
18 MONTHS

ACTUAL TREATMENT  
TIME WITH ACCELENT:  
8 MONTHS



## FINISH

Missed 2 weeks of aligner wear  
(holding pattern of night time  
only) during international trip





# CASE STUDY



## DEMOGRAPHICS

MALE  
14 YEARS OLD

## DIAGNOSIS

CLASS I, II DIV 2  
CROWDING  
DEEP BITE (50%)  
NARROW ARCHES  
EXCESSIVE OVERJET  
(5MM)

## MECHANICS

INVISALIGN®  
33 TOTAL ALIGNERS  
NO REFINEMENT  
7 DAY ALIGNER CHANGE

## TREATMENT TIME

PROJECTED TREATMENT  
TIME:  
18 MONTHS

ACTUAL TREATMENT  
TIME WITH ACCELEDENT:  
8 MONTHS



**RESULTS: 8 MONTHS**

# CASE STUDY



## DEMOGRAPHICS

MALE  
14 YEARS OLD

## DIAGNOSIS

CLASS I, II DIV 2  
CROWDING  
DEEP BITE (50%)  
NARROW ARCHES  
EXCESSIVE OVERJET (5MM)

## MECHANICS

INVISALIGN®  
33 TOTAL ALIGNERS  
NO REFINEMENT  
7 DAY ALIGNER CHANGE

## TREATMENT TIME

PROJECTED TREATMENT  
TIME:  
18 MONTHS

ACTUAL TREATMENT  
TIME WITH ACCELENT:  
8 MONTHS



## CLINICAL NOTES

Class II elastics LEFT side only to correct end-on Class II relationship and 5mm Overjet.  
Precision cutouts in the aligners at the upper left canine and lower left first molar.  
Clear button bonded directly to the upper left canine and a Caplin Hook bonded to the lower left first molar.  
Elastics were worn for the duration of treatment.

## RESULTS: 8 MONTHS

# CASE STUDY



## DEMOGRAPHICS

MALE  
14 YEARS OLD

## DIAGNOSIS

CLASS I, II DIV 2  
CROWDING  
DEEP BITE (50%)  
NARROW ARCHES  
EXCESSIVE OVERJET  
(5MM)

## MECHANICS

INVISALIGN®  
33 TOTAL ALIGNERS  
NO REFINEMENT  
7 DAY ALIGNER CHANGE

## TREATMENT TIME

PROJECTED TREATMENT  
TIME:  
18 MONTHS

ACTUAL TREATMENT  
TIME WITH ACCELENT:  
8 MONTHS



**RESULTS: 8 MONTHS**



# CASE STUDY



## DEMOGRAPHICS

MALE  
14 YEARS OLD

## DIAGNOSIS

CLASS I, II DIV 2  
CROWDING  
DEEP BITE (50%)  
NARROW ARCHES  
EXCESSIVE OVERJET  
(5MM)

## MECHANICS

INVISALIGN®  
33 TOTAL ALIGNERS  
NO REFINEMENT  
7 DAY ALIGNER CHANGE

## TREATMENT TIME

PROJECTED TREATMENT  
TIME:  
18 MONTHS

ACTUAL TREATMENT  
TIME WITH ACCELENT:  
8 MONTHS

