

CASE STUDY



DEMOGRAPHICS

FEMALE
13 YEARS OLD

DIAGNOSIS

CLASS II
SEVERE CROWDING
UNERUPTED TEETH
EXCESSIVE OVERJET

MECHANICS

INVISALGN®
44 TOTAL ALIGNERS
U: 30+14
L: 30+14
7 DAY ALIGNER CHANGE
CLASS II ELASTICS U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT TIME:
22 MONTHS

ACTUAL TREATMENT TIME WITH
ACCELEDENT:
12 MONTHS

ADDITIONAL RESOURCES

ON DEMAND CE PROGRAM

[Accelerated Aligner Therapy: Precision Taken to the Next Level](#)

Presenter: Dr. Manal Ibrahim

CASE STUDY



DEMOGRAPHICS

FEMALE
13 YEARS OLD

DIAGNOSIS

CLASS II
SEVERE CROWDING
UNERUPTED TEETH
EXCESSIVE OVERJET

MECHANICS

INVISALGN®
44 TOTAL ALIGNERS
U: 30+14
L: 30+14
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
22 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS



BASELINE



CASE STUDY



CLINCHECK® SOFTWARE TREATMENT INFORMATION

DEMOGRAPHICS

FEMALE
13 YEARS OLD

DIAGNOSIS

CLASS II
SEVERE CROWDING
UNERUPTED TEETH
EXCESSIVE OVERJET

MECHANICS


INVISALGN®
44 TOTAL ALIGNERS
U: 30+14
L: 30+14
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
22 MONTHS


ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS

TOOTH MOVEMENT ASSESMENT



Tooth	Assessment	Movement	Value	Range
UR4	Blue	Intrusion	0.5 mm	0.5 mm-1.0 mm
UR3	Blue	Root Movement	4.1 mm	4.0 mm-6.0 mm
LL4	Blue	Root Movement	4.2 mm	4.0 mm-6.0 mm
LR5	Blue	Intrusion	0.6 mm	0.5 mm-1.0 mm
Upper Right Quadrant	Black	A-P Correction	4.8 mm	> 4.0 mm
Upper Left Quadrant	Black	A-P Correction	4.2 mm	> 4.0 mm

STAGING



CASE STUDY



DEMOGRAPHICS

FEMALE
13 YEARS OLD

DIAGNOSIS

CLASS II
SEVERE CROWDING
UNERUPTED TEETH
EXCESSIVE OVERJET

MECHANICS

INVISALGN®
44 TOTAL ALIGNERS
U: 30+14
L: 30+14
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
22 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELENT:
12 MONTHS

PROGRESS Pre-Refinement



CASE STUDY



DEMOGRAPHICS

FEMALE
13 YEARS OLD

DIAGNOSIS

CLASS II
SEVERE CROWDING
UNERUPTED TEETH
EXCESSIVE OVERJET

MECHANICS

INVISALGN®
44 TOTAL ALIGNERS
U: 30+14
L: 30+14
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
22 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS



FINISH



CASE STUDY



DEMOGRAPHICS

FEMALE
13 YEARS OLD

DIAGNOSIS

CLASS II
SEVERE CROWDING
UNERUPTED TEETH
EXCESSIVE OVERJET

MECHANICS

INVISALGN®
44 TOTAL ALIGNERS
U: 30+14
L: 30+14
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
22 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS

BASELINE (2/11/16)





PRE-REFINEMENT PROGRESS 6 MONTHS (8/4/16)



FINISH 12 MONTHS (2/7/17)



CLINCHECK SOFTWARE PROJECTED OUTCOME



CASE STUDY



DEMOGRAPHICS

FEMALE
13 YEARS OLD

DIAGNOSIS

CLASS II
SEVERE CROWDING
UNERUPTED TEETH
EXCESSIVE OVERJET

MECHANICS

INVISALGN®
44 TOTAL ALIGNERS
U: 30+14
L: 30+14
7 DAY ALIGNER CHANGE
CLASS II ELASTICS
U4S/L6S

TREATMENT TIME

PROJECTED TREATMENT
TIME:
22 MONTHS

ACTUAL TREATMENT
TIME WITH ACCELEDENT:
12 MONTHS

BASELINE
(2/11/16)



PROGRESS
PRE-REFINEMENT
6 MONTHS (8/4/16)



FINISH
12 MONTHS
(2/7/17)



CLINCHECK
SOFTWARE
PROJECTED
OUTCOME

